

DNA Extraction from Cheek Cells

Cheek cell extraction

1. Scrape inside of cheek with sterile loop.
2. Twirl loop in 200 μ l of 5% Chelex buffer.
3. Add 2ul of 10mg/ml PK to tube at front bench
4. Incubate for 15-30 minutes @ 56°C
5. Vortex for 10 seconds
6. Spin at max speed for 20 seconds
7. Boil for 8 minutes in a 110°C heat block.
8. Vortex for 10 seconds
9. Spin at max speed for 2-3 minutes.
10. Use the top, aqueous solution, called the supernatant. This is where your DNA is. **STAY AWAY FROM** the pellet on the bottom where all of your cell debris and denatured proteins are. The pellet contains inhibitors that would interfere with downstream processes like quantitation or PCR.