

Alcohol Precipitation of DNA

The major considerations in using alcohols to precipitate DNA are:

Temperature: -20°C is optimal, but 0°C can be used for >20 ng/mL

2) **Amount:** For small amount of DNA (<100 ng, i.e. too little to reliably see a pellet) the use of glycogen (1 µL of a 20 mg/mL stock: Roche) will increase yield and allow visualization of the pellet.

3) **Time:** Optimal precipitation requires >1 hr at -20°C.

4) **Speed of centrifuge** (>12,000 rpm): Important for small amounts or small oligos.

5) **Mg** [10 µM] helps pellet oligos (<100 bp).

6) **Avoid precipitates:** EDTA (>10 mM) and Ca (>1 mM) will precipitate at the concentrations indicated in alcohol solutions. SDS will also precipitate with salts other than NaCl.

7) **Choice of alcohol:** Isopropanol has the advantage of requiring less volume (from 0.6 vol to 1 vol. of isopropanol is added to 1 volume of DNA/salt solution compared to 2 vol. of ethanol). However isopropanol has the disadvantage of coprecipitating more salts and is less volatile (so it air dries slower) compared to ethanol.

8) **Choice of salt:**

NaOAc, pH5.2, [0.3M]: This is the standard. (10x Stock concentration = 3 M)

NaCl [0.2M]: SDS is more soluble in NaCl (Stock concentrations vary 1M - 5M)

NH₄OAc [2 - 2.5M]: Less coprecipitation of dNTPs, good for purifying oligos. However, NH₄⁺ inhibits polynucleotide kinase).

LiCl [0.8M]: Soluble in a higher concentration of ethanol which is useful for the precipitation of RNA. However Cl⁻ inhibits the initiation of protein synthesis and RNA will display varying solubility based on size.