




AUTOMATIC MATRIX MODIFICATION for improved analytical quality and sample throughput

THGA GRAPHITE FURNACE AA



The matrix modification technique is a very important feature in the concept of interference-free trace metal determinations. Applying this technique, the chemical forms, and thereby the physical properties, of the element under study and/or the matrix, can be changed by adding a suitable reagent in excess to the sample and standard reference solutions.

Matrix modification

- Decreases the volatility of the analyte element and prevents its loss during thermal pretreatment. This also allows application of higher pretreatment temperatures for better matrix removal.
- Increases the volatility of matrix components and promotes their removal before atomization.

Perkin-Elmer furnace autosamplers feature automatic matrix modification with the ability to apply one or more matrix modifiers sequentially or simultaneously.

The following tables provide detailed information on how to prepare the most common matrix modification solutions.

Table 1. THGA Graphite Furnace: Matrix Modifiers for Routine Applications

Modifier	Major Application	Absolute Mass Required	Stock Reagents Required	Working Solutions (for a 5- μ L modifier addition)	
				Preparation	Concentration
Pd + Mg(NO ₃) ₂	Ag, As, Au, Bi, Cd, Cu, Ga, Ge, Hg, In, Mn, Sb, Se, Sn, Te, Tl	5 μ g Pd + 3 μ g Mg(NO ₃) ₂	1% (10 g/L) Pd* + 1% (10 g/L) Mg **	Dilute 1 mL of Pd stock solution and 0.1 mL of Mg stock solution to 10 mL with 18 M Ω deionized water.	0.1% (1 g/L) Pd + 0.06% (0.6 g/L) Mg(NO ₃) ₂
Pd + Mg(NO ₃) ₂	P	20 μ g Pd + 5 μ g Mg(NO ₃) ₂	1% (10 g/L) Pd* + 1% (10 g/L) Mg **	Dilute 4 mL of Pd stock solution and 0.17 mL of Mg stock solution to 10 mL with 18 M Ω deionized water.	0.4% (4 g/L) Pd + 0.1% (1 g/L) Mg(NO ₃) ₂
Mg(NO ₃) ₂	Al, Be, Co, Cr, Fe, Si	15 μ g Mg(NO ₃) ₂	1% (10 g/L) Mg*	Dilute 0.5 mL of Mg stock solution to 10 mL with 18 M Ω deionized water.	0.3% (3 g/L) Mg(NO ₃) ₂
Mg(NO ₃) ₂	Zn	5 μ g Mg(NO ₃) ₂	1% (10 g/L) Mg*	Dilute 0.17 mL of Mg stock solution to 10 mL with 18 M Ω deionized water.	0.1% (1 g/L) Mg(NO ₃) ₂

(cont'd.)

Table 1. THGA Graphite Furnace: Matrix Modifiers for Routine Applications (cont'd.)

Modifier	Major Application	Absolute Mass Required	Stock Reagents Required	Working Solutions (for a 5- μ L modifier addition)	
				Preparation	Concentration
$\text{NH}_4\text{H}_2\text{PO}_4$ + $\text{Mg}(\text{NO}_3)_2$	Cd, Pb	50 μg $\text{NH}_4\text{H}_2\text{PO}_4$ + 3 μg $\text{Mg}(\text{NO}_3)_2$	Use either: A: 10% (100 g/L) $\text{NH}_4\text{H}_2\text{PO}_4$ liquid*** + 1% (10 g/L) Mg* or B: solid ultrapure $\text{NH}_4\text{H}_2\text{PO}_4$ + 1% (10 g/L) Mg* Caution: do not use $(\text{NH}_4)_2\text{HPO}_4$	A: Dilute 1 mL of $\text{NH}_4\text{H}_2\text{PO}_4$ stock solution and 0.1 mL of Mg stock solution to 10 mL with 18 M Ω deionized water. or B: Dissolve 10 g $\text{NH}_4\text{H}_2\text{PO}_4$ in 18 M Ω deionized water and make up to 100 mL. Dilute 1 mL of this solution and 0.1 mL of Mg stock solution to 10 mL with 18 M Ω deionized water.	1.0% (10 g/L) $\text{NH}_4\text{H}_2\text{PO}_4$ + 0.06% (0.6 g/L) $\text{Mg}(\text{NO}_3)_2$

Table 2. THGA Graphite Furnace: Matrix Modifiers for Special Applications

Modifier	Major Application	Absolute Mass Required	Stock Reagents Required	Working Solutions (for a 5- μ L modifier addition)	
				Preparation	Concentration
Use either: A: 95% Ar + 5% H_2 or B: HNO_3	Halide matrices (e.g., NaCl)	Used as an alternate internal gas during dry and pyrolysis steps	None		
		NA	Ultrapure conc. HNO_3	To 80 mL 18 M Ω deionized water, add 20 mL conc. HNO_3	20% (v/v) HNO_3
CH_3OH HF HNO_3	Boric acid matrix	10 μL CH_3OH + 0.2 μL conc. HF + 0.2 μL conc. HNO_3	CH_3OH + conc. HF + conc. HNO_3	For a 10- μ L modifier addition 10 mL CH_3OH + 200 μL conc. HF + 200 μL conc. HNO_3 and mix all	

* Part No. B019-0635

** Part No. B019-0634

Note: 1% Mg corresponds to 6% $\text{Mg}(\text{NO}_3)_2$ or 10.5% $\text{Mg}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$
 1 μg $\text{Mg}(\text{NO}_3)_2$ corresponds to 0.17 μg Mg or 1.75 μg $\text{Mg}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$

*** Part No. N930-3445



PERKIN ELMER
 Look to us. And see more.

Perkin Elmer
 761 Main Avenue
 Norwalk, CT 06859-0010 USA
 Tel: 800-762-4000 or
 (+1) 203-762-4000
 Fax: (+1) 203-762-4228

Perkin Elmer is a registered trademark of affiliates of Perkin Elmer LLC.



© 1999 Perkin Elmer LLC Printed in U.S.A. ♻️

D-6124
 KG079902